



Content

- 1. General information 3
- 2. Introduction 4
- 3. Display & operating buttons 5
- 4. Preparing before operating 6
- 5. Turn on & measure 6
- 6. Calibration 8
- 7. Changing scale & temperature unit 10
- 8. Turning off 11
- 9. Cleaning & maintenance 11
- 10. Disposal 11
- 11. Technical data 12
- 12. Error codes 12
- 13. Models and scales 13

Carefully read through the operating manual even if you have prior experience with KERN refractometers.

1. General information

1.1 Intended use

The refractometer is a measuring instrument for determining the refractive index of transparent substances in liquid or in some cases also in the solid state. It is used to observe the behaviour of light as it passes from a prism with known properties to the substance being tested. Use of the refractometer for other purposes is contrary to its intended use and may be hazardous. The manufacturer shall not be liable for any damages caused by improper use.

1.2 Warranty


The warranty shall be void in the event of:
 • Failure to observe the instructions in the operating manual
 • Use for purposes other than those described
 • Modifications or opening the device housing
 • Mechanical damage and/or damage resulting from media, liquids, natural wear and tear



This digital refractometer cannot measure any liquid that is highly corrosive to metal or glass. When measuring liquids that are corrosive to plastics or react chemically with plastics, be careful not to drop the measured liquid onto the shell. Otherwise it will corrode the shell.

2. Introduction


2.1 Description



2.2 Scope of delivery
 1x Storage box | 1x Digital refractometer | 1x Operating manual | 1x AAA Battery 1.5 V | 1x Pipette | 1x Screwdriver

3. Display & operating buttons


3.1 Description display & operating buttons



Note: Please replace the battery when the [Battery icon] is displayed.


4. Preparing before operating

4.1 Install the battery



Put 1 piece of 1.5V battery into the cabin in the right way and recover the cabin again.
 Only for KERN service staff

5.1 Turn on



5.2 Measure

After turning on, clean the sample tank with distilled water and then dry it. Now fill the sample up to the mark, close the cover and press "READ".




6. Calibration

The refractometer can only be calibrated with distilled water. To do this, fill the sample tank with distilled water up to the mark and close the cover. Press "CAL" for 2 seconds to enter calibration mode. Then press "CAL" again for 2-3 seconds until "CAL" flashes in the display.

5.3 Average value measurement

Press "READ" for 2 seconds. The device starts an automatic measurement series of 15 measurements and shows the average value. Afterwards, the device automatically turns back to the normal measuring mode.




6. Calibration


Further error codes can be found in the appendix.

7. Changing scale & temperature unit

7.1 Changing scale



If the calibration was not completed successfully, an error code appears in the display. Here, for example, A01.



Press "SCALE" to change into another scales and show the converted value.

Further error codes can be found in the appendix.

7.2 Changing temperature unit



Press "SCALE" to change into another scales and show the converted value.



Press "SCALE" to change into another scales and show the converted value.

Press "SCALE" to change into another scales and show the converted value.

8. Turning off

If without any operations for 1 minute, the instrument would be automatically turned off.

9. Cleaning & maintenance

1. To avoid damages to the prism and the sample tank, clean them with distilled water after each use.
2. Dry it with a soft cloth afterwards.
3. Do not use hard or abrasive objects for cleaning.
4. Do not leave any residue in the sample tank.
5. If the refractometer is not going to be used for a longer time, remove the battery and store it at a cool and dry place.

10. Disposal

The packaging consists of environmentally friendly materials which can be disposed of via local recycling facilities. The device and storage box should be disposed

11. Technical data

Scale + accuracy + resolution	Depents to the model
Temperature	0,0 – 40,0 °C / 32,0 – 104,0 °F
Automatic Temperature Compensation	Yes
Minimum sample volume	0.2 - 0.3 ml (Marking ring)
AUTO-OFF	60 seconds
Averaging measurement	15 measurements
Battery	1 x AAA 1.5 V
Lifetime of the battery	Approx. 10.000 measurements
Overall dimensions LxWxH	125x65x30 mm
Net weight	140 g (without battery)

12. Error codes

code	Instructions
A01	Beyond the scope of calibration temperature. (0.0°C-40.0°C)
A02	During calibration, no solution or solution wrong.
A03	This instrument has a hardware failure.

13. Models and scales

Model	Scale	No.	Range	Unit	Resolution	Accuracy
ORM 10BM	Refractive Index	502	1.3330-1.4200	nD	0.0001nD	±0.0003nD
	Brix	501	0.0-99.0	%	0.1%	±0.2%
	Refractive Index	502	1.3300-1.5177	nD	0.0001nD	±0.0003nD
ORM 1RS	Refractive Index	502	1.3300-1.5177	nD	0.0001nD	±0.0003nD
	Brix	501	0.0-99.0	%	0.1%	±0.2%
	Refractive Index	502	1.3300-1.5177	nD	0.0001nD	±0.0003nD
ORM 1SU	Fructose	501	0.0-88.9	%	0.1%	±0.2%
	Glucose	502	0.0-59.9	%	0.1%	±0.2%
	Brix	503	0.0-99.0	%	0.1%	±0.2%
ORM 2SU	Refractive Index	504	1.3330-1.5177	nD	0.0001nD	±0.0003nD
	Brix	501	0.0-18.5	%	0.1%	±0.2%
	Lactose	502	0.0-15.6	%	0.1%	±0.2%
ORM 1HD	Maltose	503	0.0-19.8	%	0.1%	±0.2%
	Dextran	504	0.0-50.0	%	0.1%	±0.2%
	Brix	501	0.0-38.0	%	0.1%	±0.2%
ORM 1NA	Honey Bouma	502	0.0-99.0	%	0.1%	±0.2%
	Brix	503	0.0-99.0	%	0.1%	±0.2%
	Refractive Index	504	1.3330-1.5177	nD	0.0001nD	±0.0003nD
ORM 1SA	Salinity NaCl (%)	501	0.0-28.0	%	0.1%	±0.2%
	Salinity NaCl (°C)	502	0-28.0	°C	1%	±0.2%
	Specific Weight	503	1.000-1.220	%	0.001	±0.002
ORM 1SW	Refractive Index	504	1.3330-1.4200	nD	0.0001nD	±0.0003nD
	Salinity Seawater	501	0-100	‰	1%	±2%
	Chlorinity Seawater	502	0-37	‰	1%	±2%
ORM 1AL	Specific Weight	503	1.000-1.070	%	0.001	±0.002
	Brix	504	0.0-50.0	%	0.1%	±0.2%
	Refractive Index	505	1.3330-1.4200	nD	0.0001nD	±0.0003nD
ORM 1BR	Alcohol Meas.	501	0-72	%	1%	±1%
	Brix	503	0-50.0	%	0.1%	±0.2%
	Refractive Index	504	1.3330-1.4200	nD	0.0001nD	±0.0003nD
ORM 1W	Plato	501	0.0-30.5	°P	0.1	±0.3
	Brix	503	0.0-50.0	%	0.1%	±0.2%
	SG Wort	502	1.000-1.130	%	0.001	±0.002
ORM 2WN	Refractive Index	504	1.3330-1.4200	nD	0.0001nD	±0.0003nD
	Orchate	501	0-150	°Dm	1	±2
	Orchate France	502	0-230	°Dm	1	±2
ORM 1CD	Vol%	503	0.0-22.0	%	0.1%	±0.2%
	KMW (Baob)	501	0.0-25.0	°	0.1	±0.2
	Brix	504	0.0-50.0	%	0.1%	±0.2%
ORM 2CD	Caribe TDS 1	501	0.0-25.0	°	0.1	±0.2
	Brix	502	0.0-50.0	%	0.1%	±0.2%
	Refractive Index	503	1.3330-1.4200	nD	0.0001nD	±0.0003nD
ORM 1UN	Caribe TDS 2	501	0.00-25.00	°	0.01	±0.20
	Brix	502	0.00-30.00	%	0.01%	±0.20%
	Refractive Index	503	1.3330-1.4200	nD	0.0001nD	±0.0003nD
ORM 2UN	Urine Human	501	1.000-1.050	%	0.001	±0.002
	Urine Dog	502	1.000-1.060	%	0.001	±0.002
	Urine Cat	503	1.000-1.080	%	0.001	±0.002
ORM 1CA	Refractive Index	504	1.3330-1.4200	nD	0.0001nD	±0.0003nD
	Brix	503	0.0-50.0	%	0.1%	±0.2%
	Clear	501	1.4000-2.0	°C	0.1°C	±0.1°C
ORM 2CA	AdBlue®	502	0.0-51.0	%	0.1%	±0.2%
	Battery Fluid	503	1.000-1.050	%	0.001	±0.005
	Brix	504	0.0-50.0	%	0.1%	±0.2%
ORM 1C	Refractive Index	505	1.3330-1.4200	nD	0.0001nD	±0.0003nD
	Ethylenglycol (%)	501	0.0-100.0	%	0.1%	±0.5%
	Ethylenglycol (°C)	502	1.60-0.0	°C	0.1°C	±0.1°C
ORM 2C	Propylenglycol (%)	503	0.0-100.0	%	0.1%	±0.5%
	Propylenglycol (°C)	504	1.60-0.0	°C	0.1°C	±0.1°C
	Brix	505	0.0-90.0	%	0.1%	±0.2%