**Defender Series
Low Profile Floor Platform
Instruction Manual (DFD)**

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Table of Contents

1.	Introduction	2
1.1.	Safety Instructions.....	2
1.2.	Connecting an Indicator	3
2.	Internal structure and working principle.....	4
2.1.	Internal structure	4
2.2.	Working principle.....	5
3.	Basic requirements and installation.....	6
3.1.	Basic requirements	6
3.2.	Standard installation.....	6
3.3.	Accessory installation.....	7
4.	Operation	8
4.1.	Operation Instruction.....	8
4.2.	Precautions	8
5.	Maintenance	9
5.1.	Daily Checks	9
5.2.	Cleaning	9
5.3.	Service information	9
6.	Technical Data.....	10
6.1.	Specification	10
7.	Dimensions	11
8.	Compliance.....	12

1. Introduction

This manual contains installation, operation and maintenance instructions for the Defender Series Floor Scale.

1.1. Safety Instructions

Please follow these safety precautions:

- Position the power cord so that it does not pose a potential obstacle or disconnection hazard.
- Operate the equipment only under ambient conditions specified in the user instructions.
- Disconnect the equipment from mains power before cleaning or servicing.
- Do not operate the equipment in hazardous or explosive environments.
- The foundation at the installation site must be capable of safely supporting the weight of the floor scale plus the weight of the maximum load.
- Ensure vibrations around the scale are kept to a minimum. This also applies to transportation.
- Service should only be performed by authorized personnel.

1.2. Connecting an Indicator

When connecting the base to TD52 indicator , push the base connector onto the indicator connector and turn the locking ring clockwise to lock it in place.

When connecting the base to an indicator that does not have the matching connector, install the optional load cell cable adapter (P/N 30429931 sold separately) to the indicator. Then attach the load cell cable to the cable adapter.

Note: As an alternate wiring method, cut off the connector, strip the wires and attach them to the indicator using the wiring codes in Table 1.

The extension cable accessory (P/N: 30429931) cannot be used if the scale will be used in a legal for trade or legally controlled application, because there is no provision to seal the connection between the indicator and platform.

TABLE 1. FLOOR SCALE WIRING COLOR CODE

Indicator Signal Name	Floor Scale Wire Color
EXE+	Green
SEN+	Blue
SIG+	White
SIG-	Red
SEN-	Yellow
EXE-	Black

2. Internal structure and working principle

2.1. Internal structure


Figure2-1.Floor scale

DFD series floor scales are made up of components demonstrated above. Connect the double shielded cable to the weighing instrument and turn on the power to use. In addition, external devices such as big screens, computers and printers can be connected to the scale.

Note: The standard is single ramp.

2.2. Working principle

Put an object on the weighing platform. Under gravity, the load cell elastomer deforms and makes the impedance of the strain gauge bridge adhered to the elastomer by the strain gage out of balance. As a result, a voltage signal proportional to the weight value is outputted. The signal is then amplified by a linear amplifier, and converted to a digital signal by an analog-to-digital converter, which is processed by the microprocessor of the meter to display the weight data. Squadron staff can read the weight directly and print weight record at any time if having a printer.

The ultra-low floor platform scales have no ramp resistance and boast accurate weighing. They are suitable for use in laboratory, medical and chemical industries. Stainless steel or carbon steel can be selected according to requirements.


3. Basic requirements and installation

3.1. Basic requirements

Ensure that the installation is performed according to the instructions, and the surface is flat in the area where the scale will be positioned, The Height difference of the surface must not exceed 1mm. The maximum weight capacity of the surface at the four cell feet position should be higher than the maximum allowed weight of the cell.

3.2. Standard installation

1. Take out the weighing platform and place it on the ground.
2. Take out the footplates and place them at four corners of the platform.
3. Adjust the height of the load cell to ensure that feet are stable on the footplates.
4. Secure the ramp and mounting plate with M10 hex bolts and protective sleeves. (Check below picture for detail)


LEGEND: a. M10 hex bolts, b. Protective sleeves
Figure 3-1 Standard Installation

5. Secure four footplates with expansion bolts.
6. Connect the signal cable with the indicator.

3.3. Accessory installation

Stop bar

1. Take out the stop bar and the matching bolts.
2. On the opposite side of the ramp, use two M8 bolts to secure the stop bar to the scale body and insert rubber plugs into the two round holes.

LEGEND: c. M8 bolts, d. Rubber Plug
Figure 3-2 Optional Stop Bar Installation

Ramp

The standard is a single ramp. If you order an extra ramp, follow instruction in **section 3.2 Standard installation** to install it on the other side of the platform.


4. Operation

4.1. Operation Instruction

- After installation, if you need to use the scale for trade, you need to give the scale to local metrology department or national metrology department for verification before use.
- The scale needs to be warmed up for 30 minutes. Please refer to scale's technical manual for detailed operation.
- Metering must be performed when the "Zero Center Cursor" light on the meter is lit.

4.2. Precautions

For continued best operating performance: :

- Avoid falling loads, shocks and hard lateral impacts
- When weighing, all items should be centrally located on the floor scale and not over hanging the sides or the ramps.
- Verify the calibration of the floor scale system at regular intervals.
- Though the floor scale boast a relative strong structure, the weight placed must be within the allowed range.

5. Maintenance

5.1. Daily Checks


- Check that the scale is free from debris, especially under the edges of the scale.
- Check that there is no damage to the connecting cable.

5.2. Cleaning

Maintenance of the floor scale is limited to regular cleaning.

Clean the scale in a dry environment:

- Wipe surfaces with a damp cloth.
- Use only common household cleaning agents.
- Never use water spray or high-pressure cleaners.

When using optional Ramps or Foot Plates:

- Keep the air gap at the edge of the Floor Scale free of debris.

5.3. Service information

If error diagnosis fails to solve or describe your problem, please contact an Authorized Ohaus Service Agent. Please visit our website to locate the Ohaus office nearest you, or contact us by the information on the warranty card.

6. Technical Data

6.1. Specification

Equipment Ratings:

Indoor use only

Altitude: 2000m

Operating temperature: -10 to 40°C

Humidity: Maximum relative humidity 80% for temperatures up to 31 °C decreasing linearly to 50% relative humidity at 40°C.

Pollution degree: 2

Models	DFD600HR	DFD1500HL	DFD600WR	DFD1500WL
Capacity (kg)	600 kg	1500kg	600kg	1500kg
Maximum scale division	1:3000	1:3000	1:3000	1:3000
Junction box	Stainless steel junction box, with 6m double shielded cable			
Material	Painted steel plate, with table top thickness 6mm		304 stainless steel plate, with table top thickness 6mm	
Table top size	1000x1000mm	1250X1250mm	1000x 1000 mm	1250 x 1250 mm
Net weight	136 kg	195 kg	136 kg	195 kg

TABLE 6-1. Defender series floor scale specification

7. Dimensions



Figure 7-1. Floor scale dimensions

Models	A (mm)	B (mm)	C (mm)	D (mm)	E (mm)
DFD600HR	1000	1000	35	114	1238
DFD1500HL	1250	1250	45	114	1488
DFD600WR	1000	1000	35	114	1238
DFD1500WL	1250	1250	45	114	1488

TABLE 7-1. Floor scale dimensions

8. Compliance

Compliance to the following standards is indicated by the corresponding mark on the product.

Mark	Standard
	<p>This product complies with the applicable harmonized standards of EU Directives 2011/65/EU (RoHS). The EU Declaration of Conformity is available online.</p>
	<p>This product complies with the EU Directive 2012/19/EU (WEEE). Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment.</p> <p>For disposal instructions in Europe, refer to our website.</p>

ISO 9001 Registration

The management system governing the production of this product is ISO 9001 certified.

LIMITED WARRANTY

Ohaus products are warranted against defects in materials and workmanship from the date of delivery through the duration of the warranty period. During the warranty period Ohaus will repair, or, at its option, replace any component(s) that proves to be defective at no charge, provided that the product is returned, freight prepaid, to Ohaus.

This warranty does not apply if the product has been damaged by accident or misuse, exposed to radioactive or corrosive materials, has foreign material penetrating to the inside of the product, or as a result of service or modification by other than Ohaus. In lieu of a properly returned warranty registration card, the warranty period shall begin on the date of shipment to the authorized dealer. No other express or implied warranty is given by Ohaus Corporation. Ohaus Corporation shall not be liable for any consequential damages.

As warranty legislation differs from state to state and country to country, please contact Ohaus or your local Ohaus dealer for further details.

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