Triple Pro Balance™ Model TP2611

This manual covers Installation, Operation and Maintenance for the Ohaus Triple Pro™ Balance. To ensure proper operation of the Balance, please read this manual completely.

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Description

The Ohaus Triple Pro™ Balance is a precision weighing instrument, designed to provide years of service with virtually no maintenance. Accurate, dependable and stylish; the Ohaus Triple Pro™ offers it all, the superior quality and durability you've come to expect from Ohaus plus all the extras to satisfy every classroom activity requiring a balance. Just like the original Triple Beam from Ohaus, the Triple Pro™ gives you the confidence to know you have the accuracy and convenience of a toploading balance while allowing your students to gain hands-on experience through manipulation of the poises, attachment weights (included) and the rod & clamp assembly. The Ohaus Triple Pro™ is the perfect learning tool to help students visualize mass measurement principles and to emphasize a hands-on approach to learning.

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Unpacking and Setup

Carefully remove the balance from the carton and remove the two polystyrene inserts. Place the balance on a smooth, flat surface. Remove the separate poise (sliding weight) and three attachment weights stored in one end of the polystyrene insert and place them with your balance.

In a separate box you will find the Clamp and rod support that comes with your Triple Pro™. This is used to do specific gravity weightings. You will find a rubber washer located above the pointer. This rubber washer should be removed from the scale. Slide the separate poise up into the slot on the back of the center beam. Tilt the poise over into place on the beam. With all poises in zero position, the pointer should be near zero.

All the packing material should be retained for future use should you need to return your balance for any reason.

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Weight Storage

The three attachment weight poises should be stored in the convenient holder located beneath the beam. To do this, press in the tab on the right end of the holder and lift. The holder can now be removed. Starting at the right, place the 500g attachment weight in the slot located in the base. Place the two 1kg attachment weights into the slot in line with the 500g weight. Place the cover over the weights into the base and snap into place. Your weights are now stored for future use.

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
Zeroing Your Balance

With your balance on a smooth flat surface, with all the poises in zero position, the pointer should be near zero. For exact zero, adjust the knurled knob which is located at the left end of the beam. It is advisable to check the zero adjustment periodically during use.

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
Weighing

Place the specimen on the center of the platform and proceed as follows:

1. Starting with the largest capacity beam (500 g), move the 500 g poise to the right to the first notch which causes the pointer to drop, then, move it back one notch, causing the pointer to rise.
2. Repeat procedure with the 100 g poise.
3. Slide the 10 g poise to the position which brings the pointer to rest at zero. The weight of the specimen is the sum of the values of all poise positions, read directly from the graduated beams.

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
Introduction

Attachment Weights

Total capacity of your balance is 2610 grams when attachment weights are suspended from the pivots, and 610 grams when attachment weights are not being used.

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Care and maintenance


Keep your balance clean at all times. In general, most foreign matter may be easily removed by an air syringe, but a piece of adhesive-backed tape pressed against the magnet faces will keep them free from dirt. Never apply lubricants to knives or bearings, nor allow foreign matter to accumulate.

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Specific gravity determination:

The Model 183 Clamp and Rod Support is an accessory for elevating the balance for suspending specimens in water. The rod is inserted into the 1/2 inch (12.7 mm) recess on the underside of the base. By means of a fine wire or thread, the specimen is attached to the check pin which connects the check rod assembly to the platform loop under the base.

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Formulas:

Solids - (denser than water): To determine specific gravity of solids denser than water, the specimen should be first weighed in air and then weighed immersed in water.

Solids - (less dense than water): To determine the specific gravity of solids less dense than water, it is necessary to attach a sinker and then make the following weighings:

$$\frac{\text{Weight in air}}{\text{Weight in air minus weight in water}}$$

- a. Weight of body in air.
- b. Combined weight of body in air plus sinker in water.
- c. Weight of both body and sinker in water.

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$$\frac{a}{b - c}$$

LIMITED WARRANTY

Ohaus products are warranted against defects in, materials and workmanship from the date of delivery through the duration of the warranty period. During the warranty period Ohaus will repair, or, at its option, replace any component(s) that proves to be defective at no charge, provided that the product is returned, freight prepaid, to Ohaus.

This warranty does not apply if the product has been damaged by accident or misuse, exposed to radioactive or corrosive materials,

has foreign material penetrating to the inside of the result of service or modification by other than Ohaus. In lieu of a properly returned warranty registration card, the warranty period shall begin on the date of shipment to the authorized dealer. No other express or implied warranty is given by Ohaus Corporation. Ohaus Corporation shall not be liable for any consequential damages. As warranty legislation differs from state to state and country to country, please contact Ohaus or your local Ohaus dealer for further details.

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Printed in Poland / Impreso en Polonia

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P/N 80252244 © Ohaus Corporation 2008, all rights reserved / todos los derechos reservados.