

# IKA

designed for scientists



## **BMT-50-S Ball milling tube, stainless, 50 ml**

/// Data Sheet

Components: Tube, stainless steel balls, cover

Volume: 15 - 50 ml

Application: Dry milling of dry and brittle samples (e.g. kaolin, gypsum, coloured pigments, tablets), cell maceration, processing of materials mixed with fluids

Packing unit: 10 pieces



designed for scientists

### Technical Data

Volume range (H2O) [l]	0.02 - 0.05
Tube diameter [mm]	48
Volume Tube [ml]	50
Dry Operation	yes
FDA conformity	yes
Material: Balls	stainless steel
Material: Membrane	PP
Material: Testing Vessel and Cap	PP
Material: Gyrate Stick	PEEK
Weight [kg]	0.643
Packing unit [pcs]	10