



designed for scientists



PEN.WiCo, Wireless control button for IKA PEN electronic pipette

/// Fiche technique

The PEN.WiCo is a wireless control button that can be attached to the IKA PEN electronic pipette or used independently of the device as an alternative operating unit.

For ergonomic handling and less fatigue during repetitive tasks.



designed for scientists

Données techniques

Durée de fonctionnement admissible [%]	100
Temps de chargement [h]	1.5
Temps de fonctionnement [h]	14
Type de pile	Lithium-polymère
Capacité de la pile [mAh]	35
Distance de communication max. (en fonction de l'installation) [m]	1
Dimensions (L x H x P) [mm]	48 x 18 x 22
Poids [kg]	0.01
Plage de température du milieu admise [°C]	5 - 40
Humidité relative admissible [%]	80
Protection selon DIN EN 60529	IP 40
Interface numérique USB	USB-C
Interface WiFi	oui
Tension continue [V=]	3.7