

# IKA

designed for scientists



## **DIS-300-S-M.10 Tube de mélange avec rotor/stator, 10 pcs.**

/// Fiche technique

Volume: up to 300 ml

Application: Crushing and dissolution of tablets and dragées. This tube is equipped with a dissolver blade and opens up new application fields. It was especially designed for dissolving samples in different solvents.



designed for scientists

### Données techniques

Capacité traitée (H2O) [l]	0.1 - 0.3
Diamètre du rotor [mm]	35.5
Diamètre du tube [mm]	77
Régime de 6000 tr/min [m/s]	11.2
Volume du tube [ml]	300
Conforme à la FDA	oui
Matériau : Membrane	PP
Matériau : Tube à essai et bouchon	PP
Matériau : Tige de mélange	PEEK+PTFE+graphite+fibre de carbone
Matériau : Dissolveur	Acier inoxydable
Poids [kg]	0.66
Unité d'emballage [pcs]	10