



designed for scientists



PEN.WiCo, Wireless control button for IKA PEN electronic pipette

/// Data Sheet

The PEN.WiCo is a wireless control button that can be attached to the IKA PEN electronic pipette or used independently of the device as an alternative operating unit.

For ergonomic handling and less fatigue during repetitive tasks.



designed for scientists

Technical Data

Permissible ON time [%]	100
Charging time [h]	1.5
Working time [h]	14
Battery type	Lithium-Polymer
Battery capacity [mAh]	35
Communication distance (depend on building) max. [m]	1
Dimensions (W x H x D) [mm]	48 x 18 x 22
Weight [kg]	0.01
Permissible ambient temperature [°C]	5 - 40
Permissible relative humidity [%]	80
Protection class according to DIN EN 60529	IP 40
USB interface	USB-C
WiFi Interface	yes
DC Voltage [V=]	3.7