



designed for scientists



## DT-20 Dispersing tube, 20 ml

/// Data Sheet

Components: Tube, dispersing device, cover

Volume: 5 - 15 ml

Application: Dispersion, homogenisation, suspensions

Further suitable for:

- additional media quantities may be added during operation
- sampling during operation



designed for scientists

Packing unit: 25 pieces



designed for scientists

### Technical Data

Volume range (H2O) [l]	0.005 - 0.015
Stator diameter [mm]	22
Rotor diameter [mm]	21.5
Gap between rotor and stator [mm]	0.25
Tube diameter [mm]	29
Circumferential speed at 6000 rpm [m/s]	6.7
Volume Tube [ml]	20
FDA conformity	yes
Material: Dispersing Unit	PEEK + PTFE
Material: Membrane	PP
Material: Testing Vessel and Cap	PP
Material: Gyrate Stick	PEEK
Weight [kg]	0.432
Packing unit [pcs]	25