

IKA

designed for scientists

EN



IKA Specials Q4 / 2026
Innovative Laboratory Solutions

imLab




www.imlab.eu - info@imlab.eu



+33(0)3 20 55 19 11



+32(0)16 73 55 72



Batch mill A 10 basic

- › Grinds hard, brittle, soft and fibrous materials
- › For volumes up to 50 ml
- › Starts only when the lid is closed
- › With integrated cooling connections

Ident. No.: 0004020700

~~€ 2,380~~ € 1,666.00

-30%



2+1
Year
Warranty*



Analytical mill A 11 basic

- › Batch mill for 2 different grinding procedures
- › Impact grinding of hard, brittle or non-elastic grinding materials
- › Cutting grinding for pulverizing soft, fibrous materials with a cutting blade (not incl. with delivery)

Ident. No.: 0002900000

~~€ 2,859~~ € 2,001.30

-30%



2+1
Year
Warranty*

NEW!

Fully manufactured
by robots "Made in
Germany"



RCT basic anodized

Ident. No.: 0020133785

Magnetic Stirrer **RCT basic**

- › Heating capacity of 800 W
- › IP Protection class IP 54
- › Compact design for less space required
- › Excellent magnetic coupling at up to 1500 rpm
- › With extensive multi-level safety concept for user and application
- › Functions expandable with USB dongle (available as accessory)

Ident. No.: 0025008259

~~€ 898~~ **€ 763.30**

-15%*

*Also available as anodized version

2+8
Year
Warranty*

NEW!



RET basic anodized

Ident. No.: 0020133787

Magnetic Stirrer **RET basic**

- › Maximum speed of 1700 rpm
- › Easy connection to wireless sensors and thermometers
- › Integrated intelligent networking and digitalization features

Ident. No.: 0025008265

~~€ 1,298~~ **€ 1,103.30**

-15%*

*Also available as anodized version

2+8
Year
Warranty*

NEW!



Magnetic Stirrer **RH basic**

- › Compact, powerful magnetic stirrer for entry-level or basic laboratory use
- › External temperature control with ETS-D5 contact thermometer
- › 600 W heating power and 135 mm aluminum heating plate
- › Maximum temperature: 310 °C and maximum speed: 1500 rpm
- › All components shown are included in the scope of delivery

Ident. No.: 0020106987

~~€ 578~~ **€ 491.30**

-15%

2+1
Year
Warranty*

NEW!



Cooling and Shaking Incubator INC 125 FCS digital (SP20 | SP25)

- › High shaking speed from 10 to 300 rpm
- › Shaking stroke: 20 mm | 25 mm for optimum oxygen input
- › Decontamination mode (100 °C) for prevention of cross-contamination
- › Wide temperature range: 5 °C to 80 °C
- › Defrost function with continuous water drainage through the drain opening

Ident. No.: 0020135425 | 0020136638

~~€ 7,799~~ **€ 7,019.10**

-10%

2+1
Year
Warranty*

ECO
energy-efficient
technology



Rotary Evaporator RV 3 V-C

- › Incl. heating bath with 4 l volume for up to 99 °C
- › With single-handed manual lift for precise positioning of the glassware
- › Speed range: 20 – 300 rpm
- › High-efficiency condenser with 1500 cm² cooling surface

Ident. No.: 0010003328

~~€ 3,231~~ **€ 2,907.90**

-10%

2+3
Year
Warranty*

Compact &
user-friendly



Vacuum Technology VACSTAR lite Starter Solution

- › Includes VACSTAR lite pump + VC 10 lite controller
- › Integrated device communication with circulators for optimized workflows
- › Valve cleaning function extends system life & reduces maintenance

Ident. No.: 0010014496

~~€ 2,459~~ **€ 2,089.73**

-15%

2+1
Year
Warranty*

ECO
energy-efficient
technology

Suitable for in-vitro diagnostics



Single-Use Crushing System UTDD control Tissue Solution

- › With ULTRA-TURRAX® Tube Drive P control disperser
- › Especially designed for homogeneous crushing of tissue, cartilage or bones
- › Operated in combination with patented disposable plastic tubes
- › Closed tube with hermetical sealing to prevent cross-contamination during transport

Ident. No.: 0010013008

~~€ 1,973~~ € 1,677.05

-15%

2+1
Year
Warranty*



Overhead Stirrer RW 20 digital

- › With digital display
- › Robust, slimline, ergonomic design
- › With constant power-drive
- › Two speed ranges for universal use from 60 – 2000 rpm
- › Push-through agitator shafts (only when stationary)

Ident. No.: 0003593000

~~€ 988~~ € 889.20

-10%

2+1
Year
Warranty*



Centrifuge IKA mini G

- › Mini-centrifuge for all applications which do not require high speeds
- › Can be used for both PCR-vessels and PCR-strips
- › The centrifuge will run only when the lid is closed
- › Lid opens easily with the press of a button
- › With a fixed speed of 6000 rpm

Ident. No.: 0003958000

~~€ 378~~ € 321.30

-15%

2+1
Year
Warranty*



NEW!

Electronic Pipette IKA PEN 5 – 100 µl

2+1
Year
Warranty*

With wireless control
button PEN.WiCO
(available as accessory)

- › Available in three different volume ranges
- › Powered by a rechargeable lithium battery for operating a full working day on a single charge
- › With smart modes like reverse pipetting, multi-dispense or dilution of two liquids
- › With easy and straightforward self-calibration feature
- › Intuitive menu structure
- › Ergonomic design that fits naturally in your hand, reducing fatigue during repetitive tasks

Ident. No.: 0020025102

~~€ 399~~ **€ 359.10**

-10%



IKA PEN 0.5 – 10 µl

Ident. No.: 0020025101

IKA PEN 10 – 200 µl

Ident. No.: 0020025103

Original size

IKA PEN App

