

IKA

designed for scientists



S 25 N - 8 G - ST Outils de dispersion

/// Fiche technique

ST tools are recommended for use with tissue and other fibrous materials

- Saw tooth design creates a cutting and tearing action for shredding fibrous materials
- Made of high quality stainless steel



designed for scientists

Données techniques

Capacité traitée (H2O) [l]	0.001 - 0.05
Diamètre du sator [mm]	8
Diamètre du rotor [mm]	6.1
Interstice stator-rotor [mm]	0.25
Vitesse admissible max. [rpm]	25000
Vitesse périphérique [m/s]	8.0
Profondeur d'immersion [mm]	27 - 85
Longueur totale [mm]	108
Matériel en contact avec le produit	PTFE, AISI 316L
pH	2 - 13
Résiste aux solvants	oui
Résiste aux abrasifs	oui
Température de travail max. [°C]	180
Méthodes de stérilisation	toutes méthodes
Finesse finale, suspension [μm]	10 - 50
Finesse finale, émulsions [μm]	1 - 10
Poids [kg]	0.12