

IKA

designed for scientists



ST-50-M-gamma Stirring tube, sterile, 50 ml

/// Data Sheet

Components: Tube, stirring device, cover

Volume: 15 - 50 ml

Application: Mixing, stirring, extractions, preparation of soil sample suspensions, solubility testing for medication

Further suitable for:

- additional media quantities may be added during operation
- sampling during operation
- sample preparation under sterile conditions



- germ free conservation of samples (tissue, blood, etc.)

Packing unit: 10 pieces



designed for scientists

Technical Data

Volume range (H2O) [l]	0.015 - 0.05
Tube diameter [mm]	48
Volume Tube [ml]	50
Dry Operation	yes
FDA conformity	yes
Material: Membrane	PP
Material: Testing Vessel and Cap	PP
Material: Gyrate Stick	PEEK
Sterile	yes
Weight [kg]	0.4307
Packing unit [pcs]	10