


IKA

designed for scientists

VAN-SP-4

/// Spindle dimensions (mm)



/// Full-scale range based on speed: 1... 200 rpm

Spindle	lo-vi		me-vi		hi-vi I		hi-vi II	
	<i>M%</i>		<i>M%</i>		<i>M%</i>		<i>M%</i>	
	<i>100</i>	<i>100</i>	<i>100</i>	<i>100</i>	<i>100</i>	<i>100</i>	<i>100</i>	<i>100</i>
	Min.	Max.	Min.	Max.	Min.	Max.	Min.	Max.
VAN-SP-4	2544.8	508954	27150	5430000	54300	10860000	217200	43440000
<i>rpm</i>	<i>200</i>	<i>1</i>	<i>200</i>	<i>1</i>	<i>200</i>	<i>1</i>	<i>200</i>	<i>1</i>

/// Measuring range base on speed: 0.1... 10 rpm

Spindle	lo-vi		me-vi		hi-vi I		hi-vi II	
	<i>M%</i>		<i>M%</i>		<i>M%</i>		<i>M%</i>	
	<i>10</i>	<i>100</i>	<i>10</i>	<i>100</i>	<i>10</i>	<i>100</i>	<i>10</i>	<i>100</i>
	Min.	Max.	Min.	Max.	Min.	Max.	Min.	Max.
VAN-SP-4	5089.5	5089539	54300	54300000	108600	108600000	434400	434400000
<i>rpm</i>	<i>10</i>	<i>0.1</i>	<i>10</i>	<i>0.1</i>	<i>10</i>	<i>0.1</i>	<i>10</i>	<i>0.1</i>