



NEW



Needle Printer KERN YKT-01

High-quality printer with large LCD display and numerous functions

Technical data

- Printing process: 8 print magnets with dot matrix
- Speed: 1,7 line/s
- Character size (ASCII characters): W×H 1.7×2.6 mm
- Service life: 1 million lines
- Ribbon cartridge: EPSON ERC-22 or ERC-09B
- Operation via mains adapter 230 V/50 Hz (standard EU), as standard

Price printer without interface cable, KERN YKT-01,

Accessories

- Colour ribbon for printer (1 pc.), KERN YKT-A03,
- Paper roll (1 roll), KERN 911-013-002,
- Paper roll (10 rolls), KERN 911-013-010,

STANDARD




Easy Touch

Tip:

EasyTouch Individual Print

The print product (e.g. labels) can be individually designed. The size of the print area as well as the placement of the print elements can be customized. The barcode function enables the printing of each element as a barcode or QR code.

For further details, see *EasyTouch*



Example printouts with ET Individual Print

Important note for printer orders:

Please state your balance model when ordering, so that the printer can be preconfigured in the factory to suit your balance.

Interface Cable for RS-232

For direct connection to a printer or PC, cable length approx. 1.5 m.
The interface parameters can be set using the balance keypad.


Tip:

The interface cable is **free of charge** if you purchase balance, printer or Ethernet-/WiFi-/Bluetooth adapter together.

Printer, ...		PC, YKI-02*	911-013*, YKN-01*	YKG-01, YKH-01	YKC-01	YKT-01
Balances		KERN	KERN	KERN	KERN	KERN
440, ALS/ALJ, DE-D, EMB-V, FKA, PCD, PLS/PLJ, KDP	CDS, CKE, DS, FKB-A (until year of prod. 2022)	572-926				
EOC, IOC, BID, KIB-T, KIB-TM, UID		EOC-A12				
KGP		572-926	-			
KFB-TM, KFN-TM (until year of construction 2019)	PFB, SFB	CFS-A01				
CPB, EWJ, KFS-TM, PFB		CFS-A01				
CFS, CCS, KFS-TM, IFS (from year of production 2015)		-	CFS-A01	-	CFS-A01	
ABT-NM, PBS/PBJ		770-926 (& ABT ▶ ACS-A01)				
DBS, DLB		572-926				
DAB		MLB-A05				
ABP		MLB-A05				
ADB/ADJ, ACS/ACJ, PNS/PNJ, AXS/AXE		MLB-A05				
EG/EW, PES/PEJ		474-926 (PES/PEJ ▶ MLB-A05)				
PES/PEJ						
CCA		CCA-A01				
RPB (not compatible with YKN-01)		-	CFS-A01			
FES/FEJ, PWS		MLB-A05 (PWS-A02)				
ADS/ADT, PDS/PDT		572-926				

* ONLY WHILE STOCKS LAST

With products from KERN you are ready to face the challenges of Industry 4.0



We will help to combine industrial production with the latest information and communication technology, with the aim of improving quality, reducing costs and using time and resources more efficiently so that you are able to react to future requirements in a more flexible manner. You will benefit from our uniform data protocol combined with our BalanceConnection 4.0 data software, for details please see page 167/168.