

LED ring illumination OZB-A4571



Function

1. Connection power supply unit
2. Dimmer light intensity
3. Fixing screws



Advantages at a glance

- ✓ Constant colour temperature
- ✓ Long service life
- ✓ Maintenance-free, no lamp replacement
- ✓ Very low power consumption
Available with UK power plug



Safety informations

The light from the ring illuminator is very bright. Therefore, always switch on the illumination before looking through the eyepieces. Avoid looking directly into the LEDs.



Technical details

Illuminance : 3W-LED
 Colour temperature: 7000-11000K
 Inner Ø : 60mm
 Input voltage: AC 100-240V
 50/60HZ 0,2A max.
 Output voltage: DC 12V 0,5A



Note

The ring lighting does not fit the following series:

- OSE-1
- OSF-4G
- OZC-5
- OZG-4
- OZL-45R