

LED ring illumination OBB-A6102



Function

1. Dimmer light intensity
2. On/Off button
3. Fixing screws



Advantages at a glance

- ✓ Constant colour temperature
 - ✓ Long service life
 - ✓ Maintenance-free, no lamp replacement
 - ✓ Very low power consumption
- Available with UK power plug



Safety informations

The light from the ring illuminator is very bright. Therefore, always switch on the illumination before looking through the eyepieces. Avoid looking directly into the LEDs.



Technical details

Illuminance : 4,5W-LED
 Colour temperature: ca. 7600K
 Inner Ø: 63 mm
 Input voltage: AC 100-240V



Note

The ring lighting does not fit the following series:

OSE-1
 OSF-4G
 OZC-5
 OZG-4
 OZL-45R