



Balances & Scales Catalog



Ingenious Portfolio

Our rich history and extensive industry knowledge has resulted in a comprehensive line of high-performance balances and scales that continually meet our customers' needs at economical and cost-effective prices.

Look inside for NEW product information!

Ingeniously Practical

imLab



www.imlab.eu - info@imlab.eu



+33(0)3 20 55 19 11 +32(0)16 73 55 72

OIML E1 WEIGHTS

Wire shape, stainless steel

Model	Item N°
Weight 1 mg, OIML E1 Wire	30517704
Weight 2 mg, OIML E1 Wire	30517705
Weight 5 mg, OIML E1 Wire	30517706
Weight 10 mg, OIML E1 Wire	30517707
Weight 20 mg, OIML E1 Wire	30517958
Weight 50 mg, OIML E1 Wire	30517959
Weight 100 mg, OIML E1 Wire	30517960
Weight 200 mg, OIML E1 Wire	30517961
Weight 500 mg, OIML E1 Wire	30517962

With NVLAP Certificate

Item N°
30517963
30517964
30517965
30517966
30517967
30517968
30517969
30517970
30517971



OIML E1 WEIGHTS

Cylindrical shape, Troemner Alloy 8

Model	Item N°
Weight 1 g, OIML E1	30517972
Weight 2 g, OIML E1	30517973
Weight 5 g, OIML E1	30517974
Weight 10 g, OIML E1	30517975
Weight 20 g, OIML E1	30517976
Weight 50 g, OIML E1	30517977
Weight 100 g, OIML E1	30517978
Weight 200 g, OIML E1	30517979
Weight 500 g, OIML E1	30517980
Weight 1 kg, OIML E1	30517981
Weight 2 kg, OIML E1	30517982
Weight 5 kg, OIML E1	30517983
Weight 10 kg, OIML E1	30517984
Weight 20 kg, OIML E1	30517985

With NVLAP Certificate

Item N°
30390941
30390940
30390939
30390938
30390937
30390936
30390935
30390934
30390933
30390932
30390931
30390930
30390929
30390928



OIML E1 WEIGHTS SETS

Flat polygonal sheet, Troemner Alloy 8, stainless steel, aluminum

Model	Item N°
Weight Set (12), 500 mg-1 mg, OIML E1	30517987
Weight Set (20), 50 g-1 mg, OIML E1	30517988
Weight Set (21), 100 g-1 mg, OIML E1	30517989
Weight Set (23), 2x200 g-1 mg, OIML E1	30517990
Weight Set (24), 500 g-1 mg, OIML E1	30517991
Weight Set (25), 1 kg-1 mg, OIML E1	30517992
Weight Set (27), 2x2 kg-1 mg, OIML E1	30517993
Weight Set (28), 5 kg-1 mg, OIML E1	30517994
Weight Set (29), 10 kg-1 mg, OIML E1	30517995
Weight Set (8), 50 g-1 g, OIML E1	30517996
Weight Set (9), 100 g-1 g, OIML E1	30517997
Weight Set (11), 2x200 g-1 g, OIML E1	30517998
Weight Set (12), 500 g-1 g, OIML E1	30517999
Weight Set (13), 1 kg-1 g, OIML E1	30518000
Weight Set (15), 2x2 kg-1 g, OIML E1	30518001
Weight Set (16), 5 kg-1 g, OIML E1	30518002
Weight Set (17), 10 kg-1 g, OIML E1	30518003

With NVLAP Certificate

Item N°
30390926
30390924
30390922
30390040
30390918
30390916
30390038
30390912
30390909
30390923
30390921
30390039
30390917
30390915
30390037
30390911
30390908



All weights come with a plastic case.

Weights and weight sets with NVLAP certificates: In order to get a personalized certificate please contact our Sales Support Center on receipt of your order and provide the calibration certificate number

WEIGHT SET COMPOSITION - CLASS E1

Item N°	NVLAP Certificate	1 mg	2 mg*	5 mg	10 mg	20 mg*	50 mg	100 mg	200 mg*	500 mg	1 g	2 g*	5 g	10 g	20 g*	50 g	100 g	200 g*	500 g	1 kg	2 kg*	5 kg	10 kg
30517987		.	.	.	.	.	.	.	.	.													
30517988		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.							
30517989		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.						
30517990		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.					
30517991		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.				
30517992		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.		
30517993		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	
30517994		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
30517995		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
30517996											.	.	.	.	.	.	.	.	.	.	.	.	.
30517997											.	.	.	.	.	.	.	.	.	.	.	.	.
30517998											.	.	.	.	.	.	.	.	.	.	.	.	.
30517999											.	.	.	.	.	.	.	.	.	.	.	.	.
30518000											.	.	.	.	.	.	.	.	.	.	.	.	.
30518001											.	.	.	.	.	.	.	.	.	.	.	.	.
30518002											.	.	.	.	.	.	.	.	.	.	.	.	.
30518003											.	.	.	.	.	.	.	.	.	.	.	.	.
30390926	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
30390924	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
30390922	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
30390040	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
30390918	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
30390916	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
30390038	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
30390912	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
30390909	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
30390923	.										.	.	.	.	.	.	.	.	.	.	.	.	.
30390921	.										.	.	.	.	.	.	.	.	.	.	.	.	.
30390039	.										.	.	.	.	.	.	.	.	.	.	.	.	.
30390917	.										.	.	.	.	.	.	.	.	.	.	.	.	.
30390915	.										.	.	.	.	.	.	.	.	.	.	.	.	.
30390037	.										.	.	.	.	.	.	.	.	.	.	.	.	.
30390911	.										.	.	.	.	.	.	.	.	.	.	.	.	.
30390908	.										.	.	.	.	.	.	.	.	.	.	.	.	.

*Weight set contains two weights of this value

All weights come with a master carrying case.

- Sets where the largest weight is under 1 kg come with a brush and forceps.
- Sets where the largest weight is 1 kg – 5 kg come with brush, lifters and forceps.
- Weight sets where the largest weight is 10 kg and above are supplied with a Complete Accessory Kit which includes a microfibre cloth, nylon glove, small and wide brushes, and forceps.
- Weight and weight sets with NVLAP certificates: In order to get personalized certificate, please contact Sales Support Center.

The National Voluntary Laboratory Accreditation Program (NVLAP) provides third-party accreditation to testing and calibration laboratories in response to legislative actions or requests from government agencies or private-sector organizations. NVLAP-accredited laboratories are assessed against the management and technical requirements published in the International Standard, ISO/IEC 17025:2017. NVLAP has obtained approval from the International Laboratory Accreditation Cooperation (ILAC) for the use of the ILAC MRA Mark. For further details as well as list of member organizations please visit <https://ilac.org/>.

OIML E2 WEIGHTS

Flat polygonal sheet, aluminium, stainless steel

Model	Item N°
Weight 1 mg, OIML E2	30518004
Weight 2 mg, OIML E2	30518005
Weight 5 mg, OIML E2	30518006
Weight 10 mg, OIML E2	30518007
Weight 20 mg, OIML E2	30518008
Weight 50 mg, OIML E2	30518009
Weight 100 mg, OIML E2	30518010
Weight 200 mg, OIML E2	30518011
Weight 500 mg, OIML E2	30518012

With NVLAP Certificate

Item N°
80780284
80780285
80780286
80780287
80780288
80780289
80780290
80780291
80780292



OIML E2 WEIGHTS

Cylindrical shape, Troemner Alloy 8

Model	Item N°
Weight 1 g, OIML E2	30518013
Weight 2 g, OIML E2	30518014
Weight 5 g, OIML E2	30518015
Weight 10 g, OIML E2	30518016
Weight 20 g, OIML E2	30518017
Weight 50 g, OIML E2	30518018
Weight 100 g, OIML E2	30518019
Weight 200 g, OIML E2	30518020
Weight 500 g, OIML E2	30518021
Weight 1 kg, OIML E2	30518022
Weight 2 kg, OIML E2	30518023
Weight 5 kg, OIML E2	30518024
Weight 10 kg, OIML E2	30518025
Weight 20 kg, OIML E2	30518026

With NVLAP Certificate

Item N°
80780293
80780294
80780295
80780296
80780297
80780298
80780299
80780300
80780301
80780302
80780303
80780304
80780305
80780306



OIML E2 WEIGHTS SETS

Cylindrical shape, flat polygonal sheet Troemner Alloy 8, stainless steel, aluminum

Model	Item N°
Weight Set (12), 500 mg-1 mg, OIML E2	30518028
Weight Set (20), 50 g-1 mg, OIML E2	30518029
Weight Set (21), 100 g-1 mg, OIML E2	30518030
Weight Set (23), 2x200 g-1 mg, OIML E2	30518031
Weight Set (24), 500 g-1 mg, OIML E2	30518032
Weight Set (25), 1 kg-1 mg, OIML E2	30518033
Weight Set (27), 2x2 kg-1 mg, OIML E2	30518034
Weight Set (28), 5 kg-1 mg, OIML E2	30518035
Weight Set (8), 50 g-1 g, OIML E2	30518036
Weight Set (9), 100 g-1 g, OIML E2	30518037
Weight Set (11), 2x200 g-1 g, OIML E2	30518038
Weight Set (12), 500 g-1 g, OIML E2	30518039
Weight Set (13), 1 kg-1 g, OIML E2	30518040
Weight Set (15), 2x2 kg-1 g, OIML E2	30518041
Weight Set (16), 5 kg-1 g, OIML E2	30518042

With NVLAP Certificate

Item N°
80780381
80780383
80780385
80780387
80780391
80780393
80780395
80780399
80780384
80780386
80780388
80780392
80780394
80780396
80780400



All weights come with a plastic case.

Weights and weight sets with NVLAP certificates: In order to get a personalized certificate please contact our Sales Support Center on receipt of your order and provide the calibration certificate number

WEIGHT SET COMPOSITION - CLASS E2

Item N°	NVLAP Certificate	1 mg	2 mg*	5 mg	10 mg	20 mg*	50 mg	100 mg	200 mg*	500 mg	1 g	2 g*	5 g	10 g	20 g*	50 g	100 g	200 g*	500 g	1 kg	2 kg*	5 kg	10 kg
30518028		.	.	.	.	.	.	.	.	.													
30518029		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.							
30518030		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.						
30518031		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.					
30518032		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.				
30518033		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.		
30518034		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	
30518035		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
30518036											.	.	.	.	.	.	.	.	.	.	.	.	.
30518037											.	.	.	.	.	.	.	.	.	.	.	.	.
30518038											.	.	.	.	.	.	.	.	.	.	.	.	.
30518039											.	.	.	.	.	.	.	.	.	.	.	.	.
30518040											.	.	.	.	.	.	.	.	.	.	.	.	.
30518041											.	.	.	.	.	.	.	.	.	.	.	.	.
30518042											.	.	.	.	.	.	.	.	.	.	.	.	.
80780381	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
80780383	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
80780385	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
80780387	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
80780391	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
80780393	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
80780395	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
80780399	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
80780384	.										.	.	.	.	.	.	.	.	.	.	.	.	.
80780386	.										.	.	.	.	.	.	.	.	.	.	.	.	.
80780388	.										.	.	.	.	.	.	.	.	.	.	.	.	.
80780392	.										.	.	.	.	.	.	.	.	.	.	.	.	.
80780394	.										.	.	.	.	.	.	.	.	.	.	.	.	.
80780396	.										.	.	.	.	.	.	.	.	.	.	.	.	.
80780400	.										.	.	.	.	.	.	.	.	.	.	.	.	.

*Weight set contains two weights of this value

All weights come with a master carrying case.

- Sets where the largest weight is under 1 kg come with a brush and forceps.
- Sets where the largest weight is 1 kg – 5 kg come with brush, lifters and forceps.
- Weight sets where the largest weight is 10 kg and above are supplied with a Complete Accessory Kit which includes a microfibre cloth, nylon glove, small and wide brushes, and forceps.
- Weight and weight sets with NVLAP certificates: In order to get personalized certificate, please contact Sales Support Center.

The National Voluntary Laboratory Accreditation Program (NVLAP) provides third-party accreditation to testing and calibration laboratories in response to legislative actions or requests from government agencies or private-sector organizations. NVLAP-accredited laboratories are assessed against the management and technical requirements published in the International Standard, ISO/IEC 17025:2017. NVLAP has obtained approval from the International Laboratory Accreditation Cooperation (ILAC) for the use of the ILAC MRA Mark. For further details as well as list of member organizations please visit <https://ilac.org/>.

OIML F1 WEIGHTS

Flat polygonal sheet, aluminium, stainless steel

Model	Item N°
Weight, 1 mg, OIML F1	80780307
Weight, 2 mg, OIML F1	80780308
Weight, 5 mg, OIML F1	80780309
Weight, 10 mg, OIML F1	80780310
Weight, 20 mg, OIML F1	80780311
Weight, 50 mg, OIML F1	80780312
Weight, 100 mg, OIML F1	80780313
Weight, 200 mg, OIML F1	80780314
Weight, 500 mg, OIML F1	80780315

With NVLAP Certificate

Item N°
80780331
80780332
80780333
80780334
80780335
80780336
80780337
80780338
80780339



OIML F1 WEIGHTS

Cylindrical shape, Troemner Alloy 8, stainless steel

Model	Item N°
Weight, 1 g, OIML F1	80780316
Weight, 2 g, OIML F1	80780317
Weight, 5 g, OIML F1	80780318
Weight, 10 g, OIML F1	80780319
Weight, 20 g, OIML F1	80780320
Weight, 50 g, OIML F1	80780321
Weight, 100 g, OIML F1	80780322
Weight, 200 g, OIML F1	80780323
Weight, 500 g, OIML F1	80780324
Weight, 1 kg, OIML F1	80780325
Weight, 2 kg, OIML F1	80780326
Weight, 5 kg, OIML F1	80780327
Weight, 10 kg, OIML F1	80780328
Weight, 20 kg, OIML F1	80780329

With NVLAP Certificate

Item N°
80780340
80780341
80780342
80780343
80780344
80780345
80780346
80780347
80780348
80780349
80780350
80780351
80780352
80780353



OIML F1 WEIGHTS SETS

Cylindrical shape, flat polygonal sheet Troemner Alloy 8, stainless steel, aluminum

Model	Item N°
Weight Set (12), 500 mg-1 mg, OIML F1	80780408
Weight Set (20), 50 g-1 mg, OIML F1	80780410
Weight Set (21), 100 g-1 mg, OIML F1	80780412
Weight Set (23), 2x200 g-1 mg, OIML F1	80780414
Weight Set (24), 500 g-1 mg, OIML F1	80780418
Weight Set (25), 1 kg-1 mg, OIML F1	80780420
Weight Set (27), 2x2 kg-1 mg, OIML F1	80780422
Weight Set (28), 5 kg-1 mg, OIML F1	80780426
Weight Set (8), 50 g-1 g, OIML F1	80780411
Weight Set (9), 100 g-1 g, OIML F1	80780413
Weight Set (11), 2x200 g-1 g, OIML F1	80780415
Weight Set (12), 500 g-1 g, OIML F1	80780419
Weight Set (13), 1 kg-1 g, OIML F1	80780421
Weight Set (15), 2x2 kg-1 g, OIML F1	80780423
Weight Set (16), 5 kg-1 g, OIML F1	80780427

With NVLAP Certificate

Item N°
80780435
80780437
80780439
80780441
80780445
80780447
80780449
80780453
80780438
80780440
80780442
80780446
80780448
80780450
80780454



All weights come with a plastic case.

Weights and weight sets with NVLAP certificates: In order to get a personalized certificate please contact our Sales Support Center on receipt of your order and provide the calibration certificate number

WEIGHT SET COMPOSITION - CLASS F1

Item N°	NVLAP Certificate	1 mg	2 mg*	5 mg	10 mg	20 mg*	50 mg	100 mg	200 mg*	500 mg	1 g	2 g*	5 g	10 g	20 g*	50 g	100 g	200 g*	500 g	1 kg	2 kg*	5 kg	10 kg
80780408		.	.	.	.	.	.	.	.	.													
80780410		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.							
80780412		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.						
80780414		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.					
80780418		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.				
80780420		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.		
80780422		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	
80780426		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
80780411											.	.	.	.	.	.	.	.	.	.	.	.	.
80780413											.	.	.	.	.	.	.	.	.	.	.	.	.
80780415											.	.	.	.	.	.	.	.	.	.	.	.	.
80780419											.	.	.	.	.	.	.	.	.	.	.	.	.
80780421											.	.	.	.	.	.	.	.	.	.	.	.	.
80780423											.	.	.	.	.	.	.	.	.	.	.	.	.
80780427											.	.	.	.	.	.	.	.	.	.	.	.	.
80780435	.	.	.	.	.	.	.	.	.	.													
80780437	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
80780439	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
80780441	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
80780445	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
80780447	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
80780449	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
80780453	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
80780438	.										.	.	.	.	.	.	.	.	.	.	.	.	.
80780440	.										.	.	.	.	.	.	.	.	.	.	.	.	.
80780442	.										.	.	.	.	.	.	.	.	.	.	.	.	.
80780446	.										.	.	.	.	.	.	.	.	.	.	.	.	.
80780448	.										.	.	.	.	.	.	.	.	.	.	.	.	.
80780450	.										.	.	.	.	.	.	.	.	.	.	.	.	.
80780454	.										.	.	.	.	.	.	.	.	.	.	.	.	.

*Weight set contains two weights of this value.

All weights come with a master carrying case.

- Sets where the largest weight is under 1 kg come with a brush and forceps.
- Sets where the largest weight is 1 kg – 5 kg come with brush, lifters and forceps.
- Weight sets where the largest weight is 10 kg and above are supplied with a Complete Accessory Kit which includes a microfibre cloth, nylon glove, small and wide brushes, and forceps.
- Weight and weight sets with NVLAP certificates: In order to get personalized certificate, please contact Sales Support Center.

The National Voluntary Laboratory Accreditation Program (NVLAP) provides third-party accreditation to testing and calibration laboratories in response to legislative actions or requests from government agencies or private-sector organizations. NVLAP-accredited laboratories are assessed against the management and technical requirements published in the International Standard, ISO/IEC 17025:2017. NVLAP has obtained approval from the International Laboratory Accreditation Cooperation (ILAC) for the use of the ILAC MRA Mark. For further details as well as list of member organizations please visit <https://ilac.org/>.

OIML F2 WEIGHTS

Flat polygonal sheet, aluminium, stainless steel

Model	Item N°
Weight, 1 mg, OIML F2	30390876
Weight, 2 mg, OIML F2	30390875
Weight, 5 mg, OIML F2	30390874
Weight, 10 mg, OIML F2	30390873
Weight, 20 mg, OIML F2	30390872
Weight, 50 mg, OIML F2	30390871
Weight, 100 mg, OIML F2	30390870
Weight, 200 mg, OIML F2	30390869
Weight, 500 mg, OIML F2	30390868

With NVLAP Certificate

Item N°
30390900
30390899
30390898
30390897
30390896
30390895
30390894
30390893
30390892



OIML F2 WEIGHTS

Cylindrical shape, Troemner Alloy 8, stainless steel

Model	Item N°
Weight, 1 g, OIML F2	30390867
Weight, 2 g, OIML F2	30390866
Weight, 5 g, OIML F2	30390865
Weight, 10 g, OIML F2	30390864
Weight, 20 g, OIML F2	30390863
Weight, 50 g, OIML F2	30390862
Weight, 100 g, OIML F2	30390861
Weight, 200 g, OIML F2	30390860
Weight, 500 g, OIML F2	30390859
Weight, 1 kg, OIML F2	30390858
Weight, 2 kg, OIML F2	30390857
Weight, 5 kg, OIML F2	30390856
Weight, 10 kg, OIML F2	30390855
Weight, 20 kg, OIML F2	30390854

With NVLAP Certificate

Item N°
30390891
30390890
30390889
30390888
30390887
30390886
30390885
30390884
30390883
30390882
30390881
30390880
30390879
30390878



OIML F2 WEIGHTS SETS

Cylindrical shape, flat polygonal sheet Troemner Alloy 8, Stainless Steel, Aluminum

Model	Item N°
Weight Set (12), 500 mg-1 mg, OIML F2	30390829
Weight Set (20), 50 g-1 mg, OIML F2	30390827
Weight Set (21), 100 g-1 mg, OIML F2	30390825
Weight Set (23), 2x200 g-1 mg, OIML F2	30390032
Weight Set (24), 500 g-1 mg, OIML F2	30390821
Weight Set (25), 1 kg-1 mg, OIML F2	30390819
Weight Set (27), 2x2 kg-1 mg, OIML F2	30390030
Weight Set (28), 5 kg-1 mg, OIML F2	30390815
Weight Set (8), 50 g-1 g, OIML F2	30390826
Weight Set (9), 100 g-1 g, OIML F2	30390824
Weight Set (11), 2x200 g-1 g, OIML F2	30390031
Weight Set (12), 500 g-1 g, OIML F2	30390820
Weight Set (13), 1 kg-1 g, OIML F2	30390818
Weight Set (15), 2x2 kg-1 g, OIML F2	30390029
Weight Set (16), 5 kg-1 g, OIML F2	30390814

With NVLAP Certificate

Item N°
30390852
30390850
30390848
30390036
30390844
30390842
30390034
30390838
30390849
30390847
30390035
30390843
30390841
30390033
30390837



All weights come with a plastic case.

Weights and weight sets with NVLAP certificates: In order to get a personalized certificate please contact our Sales Support Center on receipt of your order and provide the calibration certificate number

WEIGHT SET COMPOSITION - CLASS F2

Item N°	NVLAP Certificate	1 mg	2 mg*	5 mg	10 mg	20 mg*	50 mg	100 mg	200 mg*	500 mg	1 g	2 g*	5 g	10 g	20 g*	50 g	100 g	200 g*	500 g	1 kg	2 kg*	5 kg	10 kg
30390829		.	.	.	.	.	.	.	.	.													
30390827		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.							
30390825		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.						
30390032		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.					
30390821		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.				
30390819		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.		
30390030		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	
30390815		.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
30390826											.	.	.	.	.	.	.	.	.	.	.	.	.
30390824											.	.	.	.	.	.	.	.	.	.	.	.	.
30390031											.	.	.	.	.	.	.	.	.	.	.	.	.
30390820											.	.	.	.	.	.	.	.	.	.	.	.	.
30390818											.	.	.	.	.	.	.	.	.	.	.	.	.
30390029											.	.	.	.	.	.	.	.	.	.	.	.	.
30390814											.	.	.	.	.	.	.	.	.	.	.	.	.
30390852	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
30390850	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
30390848	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
30390036	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
30390844	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
30390842	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
30390034	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
30390838	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
30390849	.										.	.	.	.	.	.	.	.	.	.	.	.	.
30390847	.										.	.	.	.	.	.	.	.	.	.	.	.	.
30390035	.										.	.	.	.	.	.	.	.	.	.	.	.	.
30390843	.										.	.	.	.	.	.	.	.	.	.	.	.	.
30390841	.										.	.	.	.	.	.	.	.	.	.	.	.	.
30390033	.										.	.	.	.	.	.	.	.	.	.	.	.	.
30390837	.										.	.	.	.	.	.	.	.	.	.	.	.	.

*Weight set contains two weights of this value.

All weights come with a master carrying case.

- Sets where the largest weight is under 1 kg come with a brush and forceps.
- Sets where the largest weight is 1 kg – 5 kg come with brush, lifters and forceps.
- Weight sets where the largest weight is 10 kg and above are supplied with a Complete Accessory Kit which includes a microfibre cloth, nylon glove, small and wide brushes, and forceps.
- Weight and weight sets with NVLAP certificates: In order to get personalized certificate, please contact Sales Support Center.

The National Voluntary Laboratory Accreditation Program (NVLAP) provides third-party accreditation to testing and calibration laboratories in response to legislative actions or requests from government agencies or private-sector organizations. NVLAP-accredited laboratories are assessed against the management and technical requirements published in the International Standard, ISO/IEC 17025:2017. NVLAP has obtained approval from the International Laboratory Accreditation Cooperation (ILAC) for the use of the ILAC MRA Mark. For further details as well as list of member organizations please visit <https://ilac.org/>.

OIML M1 WEIGHTS

Stainless steel

Model	Item N°
Weight 1 g OIML M1	11117141
Weight 2 g OIML M1	11117142
Weight 5 g OIML M1	11117143
Weight 10 g OIML M1	11117144
Weight 20 g OIML M1	11117145
Weight 50 g OIML M1	11117146
Weight 100 g OIML M1	11117147
Weight 200 g OIML M1	11117148
Weight 500 g OIML M1	11117149
Weight 1 kg OIML M1	11117150
Weight 2 kg OIML M1	11117151
Weight 5 kg OIML M1	11117152
Weight 10 kg OIML M1	11117153

Iron

Model	Item N°
Rectangular weight 20 kg M1	87000007
Rectangular weight 10 kg M1	87000008
Rectangular weight 5 kg M1	87000009

With certificate

Model	Item N°
Weight OH 1 g M1AC	11117220
Weight OH 2 g M1AC	11117221
Weight OH 5 g M1AC	11117222
Weight OH 10 g M1AC	11117223
Weight OH 20 g M1AC	11117224
Weight OH 50 g M1AC	11117225
Weight OH 100 g M1AC	11117226
Weight OH 200 g M1AC	11117227
Weight OH 500 g M1AC	11117228
Weight OH 1 kg M1AC	11117229
Weight OH 2 kg M1AC	11117230
Weight OH 5 kg M1AC	11117231
Weight OH 10 kg M1AC	11117232

OIML M3 WEIGHTS

Stainless steel

Model	Item N°
Weight 1 g OIML M3	11117206
Weight 2 g OIML M3	11117207
Weight 5 g OIML M3	11117208
Weight 10 g OIML M3	11117209
Weight 20 g OIML M3	11117210
Weight 50 g OIML M3	11117211
Weight 100 g OIML M3	11117212
Weight 200 g OIML M3	11117213
Weight 500 g OIML M3	11117214
Weight 1 kg OIML M3	11117215
Weight 2 kg OIML M3	11117216
Weight 5 kg OIML M3	11117217
Weight 10 kg OIML M3	11117218

With certificate

Model	Item N°
Weight 1 g OIML M3	11117234
Weight 2 g OIML M3	11117235
Weight 5 g OIML M3	11117236
Weight 10 g OIML M3	11117237
Weight 20 g OIML M3	11117238
Weight 50 g OIML M3	11117239
Weight 100 g OIML M3	11117240
Weight 200 g OIML M3	11117241
Weight 500 g OIML M3	11117242
Weight 1 kg OIML M3	11117243
Weight 2 kg OIML M3	11117244
Weight 5 kg OIML M3	11117245
Weight 10 kg OIML M3	11117246

OIML Weight sets 1 mg

Model	Item N°
M1 Weight set 1 mg - 500 mg (ALU)	87000010
M1 Weight set 1 mg - 100 g (ALU/Brass, SST)	87000011
M1 Weight set 1 mg - 200 g (ALU/Brass, SST)	87000012

With certificate

Model	Item N°
F1 Weight set 1 mg – 500 mg	87000034

OIML Weight sets 1 g

Model	Item N°
M1 Weight set 1 g - 50 g (Brass)	87000013
M1 Weight set 1 g - 500 g (SST, Brass)	87000014
M1 Weight set 1 g - 1000 g (SST, Brass)	87000015
M1 Weight set 1 g - 2000 g (SST, Brass)	87000016

With certificate

Model	Item N°
F1 Weight set 1 g – 500 g	87000035