

stakpure ion exchange system DS 750

stakpure

stakpure pressureless ionexchange systems, the economic way for producing demineralized water. Inside of the cartridge water flows through a bed of especially mixed anionic- and cationic resins and so is purified by the well-proven method of ion exchanging. Complete system comes with an analogue conductivity meter, wall holder and tubing kit. Designed for doctors surgeries, small labs or for gaining battery water up to 15 liters a day.

Scope of delivery and technical data

Pressure less plastic cartridge

Filling of high-quality mixed bed resins

Analogue conductivity meter

Type	DS 750
Material	Polyethylene
Flow rate	100 l/h
Capacity (at 10°dH - total salt)	750 l
Pure water quality	0.1 – 20 µS/cm
Maximum water temperature	30°C
Electrical connection	230 V / 50-60 Hz
Diameter x height (replacement)	175 x 470 (350) mm
Connectors	R ¾" male thread
Approximate weight (replacement)	6 (5) kg



Article numbers

11170500	Demineralising system DS 750
11175000	Replacement cartridge for DS 750
13140400	Wall-mounting