

CORG51



NFX
15-211



Activated carbon filter for organic and corrosive vapors - for filtration box H50C (or old cabinets range 12.X before 2002)

COMPLIANCE

ADVANTAGES

- Specific filtration of the chemicals
- NFX 15-211 compliant
- Storage of flammables, corrosives and toxics



TECHNICAL FEATURES

Model	CORG51			HS Code	38021000		
Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Storage capacity (L)	Containment sump	Shelves	Load capacity	Weight (kg)
CORG51	70 x 390 x 390		0		0	0	9.00