

FOR50



NFX
15-211



Activated carbon filter for formaldehyde vapors - for filtration box H50C (or old cabinets range 12.X before 2002)

COMPLIANCE

ADVANTAGES

- Specific filtration of the chemicals
- NF X 15-211 compliant
- Storage of flammables, corrosives and toxics



TECHNICAL FEATURES

Model	FOR50			HS Code			
Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Storage capacity (L)	Containment sump	Shelves	Load capacity	Weight (kg)
FOR50	70 x 390 x 390		0		0	0	8.00