

H061D+

NFX  
15-211prEN  
17242

# Ductless fume hood 600 mm - 1 Module - two openings

## COMPLIANCE

NF X 15-211 Classe 2

## ADVANTAGES

- CE certified
- NF X 15-211 compliant
- Storage of flammables, corrosives and toxics
- Draft standard EN 17242

## BUT ALSO...

- NF X 15-211 compliant
- Delivered ready to use
- Lighting control allowing to survey the proper working of the fan
- Low acoustic level
- Front control window to see immediately the filter
- Glass side walls
- Ergonomic hand traps or trapezoidal opening
- Holes for the cables and fluid going through
- Easy replacement of the active coal filter



## TECHNICAL FEATURES

Model	<b>H061D+</b>	Hour counter	<b>included</b>
Compliance	<b>NF X 15-211 Classe 2</b>	Sampling outlet	<b>included</b>
Airflow (m/s)	<b>0,5</b>	Backwall transparent	<b>included</b>
Voltage (volts)	<b>220-240</b>	Melamine work-top	<b>option</b>
Flow (m <sup>3</sup> /h)	<b>170</b>	Electrical power (W)	<b>63</b>
Openings for passages of cables	<b>4</b>	Sound power of the fan dB(A)	<b>60.00</b>
Anemometer	<b>included</b>	HS Code	<b>84146000</b>
Guide of chemicals adsorbed	<b>included</b>		

## COMPLEMENTARY PRODUCTS

- Activated carbon filter for ammonia vapors - for old fume hoods type 200, 201, 210 and 211 (Réf : AMM200)
- Kit to use the fume hood without activated carbon filter - for old fume hoods type 200, 201, 210 and 211 (Réf : CAR200)
- Activated carbon filter for organic and corrosive vapors - for old fume hoods type 200, 201, 210 and 211 (Réf : CORG201)
- Activated carbon filter for formaldehyde vapors - for old fume hoods type 200, 201, 210 and 211 (Réf : FOR200)
- HEPA filter for ductless fume hoods H201 - H211 (Réf : HEP200)
- HEPA filter frame for H201 - H211 (Réf : SH200)
- Containment sump and melamine work-top for H06 (Réf : BB06)
- Containment sump and tempered glass work-top for H06 (Réf : BV06)

## OPTIONS

- Lighting for ductless fume hoods (Réf : LUMI)
- Table-cabinet for H06 (Réf : TA06)
- Trolley with adjustable shelf for H06 (Réf : TR06)
- Tubular table with feet for model H06 (Réf : TTF06)
- Tubular table with wheels for model H06 (Réf : TTR06)

Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Weight (kg)
H061D+	950 x 600 x 650	735 x 545 x 505	23.00



**imLab**  
SCIENTIFIC EQUIPMENT  
FOR LABORATORY & INDUSTRY

@ Web contact

www.imlab.eu  
info@imlab.eu

☎ Phone

+33(0)3 20 55 19 11  
+32(0)16 73 55 72