

H70



NF X 15-211 Mini ductless fume hood

COMPLIANCE

NF X 15-211

ADVANTAGES

- CE certified
- NF X 15-211 compliant
- Storage of flammables, corrosives and toxics

BUT ALSO...

- Containment sump integrated
- Easy replacement of the active coal filter
- Low acoustic level
- Front control window to see immediately the filter
- NF X 15-211 compliant
- Delivered ready to use
- Lighting control allowing to survey the proper working of the fan
- Glass side walls
- Ergonomic hand traps or trapezoidal opening
- Adjustable leveling feet for enhance stability



TECHNICAL FEATURES

Model	H70	Hour counter	Oui
Compliance	NF X 15-211	Sampling outlet	Oui
Airflow (m/s)	0,5	Backwall transparent	Non
Voltage (volts)	220-240	Electrical power (W)	63
Flow (m3/h)	180	Sound power of the fan dB(A)	60.00
Anemometer	Oui	HS Code	84146000
Guide of chemicals adsorbed	Oui		

COMPLEMENTARY PRODUCTS

- Activated carbon filter for ammonia vapors - for old fume hoods type 200, 201, 210 and 211 (Réf : AMM200)
- Activated carbon filter for organic and corrosive vapors - for filtration box H50C (or old cabinets range 12.X before 2002) (Réf : CORG51)
- Activated carbon filter for formaldehyde vapors - for filtration box H50C (or old cabinets range 12.X before 2002) (Réf : FOR50)
- Activated carbon filter for organic vapors - for filtration box H50C (or old cabinets range 12.X before 2002) (Réf : ORG50)

Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Weight (kg)
H70+S	605 x 505 x 485	410 x 455 x 445	25.00



@ Web contact

www.imlab.eu
info@imlab.eu

☎ Phone

+33(0)3 20 55 19 11
+32(0)16 73 55 72