

B235



Containment sump for 2- door cabinets Types 30, 60 & 90

COMPLIANCE

ADVANTAGES

- 30-minute fire resistant and made in France
- EN 14470-1 and FM 6050 compliant
- Storage of flammables, corrosives and toxics



TECHNICAL FEATURES

Model **B235** HS Code **94032080**
Shipping dimensions HxWxD (mm) **1000 x 500 x 500**

| Model | External dimensions HxWxD (mm) | Internal dimensions HxWxD (mm) | Storage capacity (L) | Containment sump | Shelves | Load capacity | Weight (kg) |
|-------|--------------------------------|--------------------------------|----------------------|------------------|---------|---------------|-------------|
| B235 | 100 x 1000 x 410 | | 0 | | 0 | 0 | 10.00 |