

3032E



2- door underbench safety cabinet Type 30

COMPLIANCE

EN 14470-1, FM 6050

ADVANTAGES

- 30-minute fire resistant and made in France
- EN 14470-1 and FM 6050 compliant
- Storage of flammables

BUT ALSO...

- Thermo-expanding seals of the doors
- Adjustable leveling feet for enhance stability
- Holes for the connection to the ventilation
- Normalized labelling provided with the equipment
- Pedestal area under-driveable by a forklift by removing panel
- Tested and approved according to FM 6050



TECHNICAL FEATURES

| | | | |
|-------------------|--------------------------------|--------------------------------|--|
| Model | 3032E | Shipping dimensions HxWxD (mm) | 835 x 1200 x 800 |
| Compliance | EN 14470-1, FM 6050 | Ventilation | Both inlet and outlet on the back wall of the cabinet |
| Cabinet materials | Melamine wood structure | Openings diameter (mm) | 100 |
| Doors color | light grey RAL 7035 | Number of compartments | 1 |
| Number of doors | 2 | Containment sump | 1 |
| Doors type | swing doors | Grating | 1 |
| Hinges | right & left | HS Code | 94036000 |
| Door closure | self-closing | | |

ADDITIONAL ACCESSORIES RECOMMENDED - FILTRATION OR VENTILATION BOXES

- Ventilation box for a connection to the outside (Réf : CDV-A)
- Connecting kit (Réf : KRC)
- Vertical recycling air box (delivered without filter) (Réf : H40+)

OPTIONS

- Containment sump for 2- door cabinets Types 30, 60 & 90 (Réf : B235)
- Perforated cover plate for 2-door safety cabinets Types 30, 60 & 90 (Réf : C235)
- Retention shelf for 2-door safety cabinets Types 30, 60 & 90 (Réf : E35)
- Containment sump in PVC for 1-door or multirisk safety cabinets Types 30, 60 & 90 (Réf : SB148)
- Protection shelf in PVC for 1-door or multirisk safety cabinets Types 30, 60 & 90 (Réf : SE48)

| Model | External dimensions HxWxD (mm) | Internal dimensions HxWxD (mm) | Storage capacity (L) | Containment sump | Shelves | Load capacity | Weight (kg) |
|-------|--------------------------------|--------------------------------|----------------------|------------------|---------|---------------|-------------|
| 3032E | 635 x 1115 x 630 | 365 x 1050 x 420 | 60 | 1 | 0 | 50 | 108.00 |

