

# 1- door tall safety cabinet

## COMPLIANCE

## ADVANTAGES

- Made in France
- Complies with risk prevention and safety recommendations
- Storage of flammables, corrosives and toxics

### BUT ALSO...

- Complies with risk prevention and safety recommendations
- Holes for the connection to the ventilation
- Locking system
- Normalized labelling provided with the equipment



## TECHNICAL FEATURES

Model	<b>A1T</b>	Shipping dimensions HxWxD (mm)	<b>2025 x 1200 x 800</b>
Cabinet materials	<b>10/10th epoxy painted steel structure</b>	Ventilation	<b>Inlet located at the bottom of the cabinet, outlet located at the top</b>
Doors color	<b>white RAL 9010</b>	Openings diameter (mm)	<b>100</b>
Number of doors	<b>1</b>	Number of compartments	<b>1</b>
Doors type	<b>sliding doors</b>	HS Code	<b>94032080</b>
Door closure	<b>manual closing</b>		

## ADDITIONAL ACCESSORIES RECOMMENDED - FILTRATION OR VENTILATION BOXES

- Ventilation box for a connection to the outside (Réf : CDV-A)
- Connecting kit (Réf : KRC)

## OPTIONS

- Retention shelf for A1T - A2T - A3T - A4T (Réf : ES1T)

Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Storage capacity (L)	Containment sump	Shelves	Load capacity	Weight (kg)
A1T	1825 x 560 x 545	1800 x 550 x 475	120		4	50	80.00