

4- door tall safety cabinet - 4 compartments, including PVC equipments for acids

COMPLIANCE

ADVANTAGES

- Made in France
- Complies with risk prevention and safety recommendations
- Storage of flammables, corrosives and toxics

BUT ALSO...

- PVC protections for acid compartment (included)
- Allows the separation of different types of chemicals
- Normalized labelling provided with the equipment
- Holes for the connection to the ventilation
- Complies with risk prevention and safety recommendations
- Locking system



TECHNICAL FEATURES

Model	A2+2PCA	Ventilation	Inlet located at the bottom of the cabinet, outlet located at the top 100
Cabinet materials	10/10th epoxy painted steel structure	Openings diameter (mm)	
Doors color	white RAL 9010	Number of compartments	4
Number of doors	4	Shelves in PVC	2
Doors type	swing doors	Containment sump	3
Hinges	right	Containment sump in PVC	1
Door closure	manual closing	HS Code	94032080
Shipping dimensions HxWxD (mm)	2025 x 1200 x 800		

ADDITIONAL ACCESSORIES RECOMMENDED - FILTRATION OR VENTILATION BOXES

- Ventilation box for a connection to the outside (Réf: CDV-A)
- Connecting kit (Réf: KRC)

OPTIONS

- Retention shelf (Réf: EE1)
- Containment sump in PVC (Réf: SBPVC)
- Shelf protection in PVC (Réf: SEPVC)

Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Storage capacity (L)	Containment sump	Shelves	Load capacity	Weight (kg)
A2+2PCA	1825 x 1100 x 545	4x (890 x 550 x 475)	220	3	8	50	170.00