

APC112+



# 2- door working cover melamine safety cabinet for acids and bases

## COMPLIANCE

## ADVANTAGES

- Made in France
- Complies with risk prevention and safety recommendations
- Storage of corrosives

### BUT ALSO...

- Adjustable leveling feet for enhance stability
- Adjustable leveling feet for enhance stability
- Cabinet with compartments to separate acids and bases
- PP containment sumps on sliding rails
- Complies with risk prevention and safety recommendations
- Built in melamine with a great resistance to corrosive vapors
- Holes for the connection to the ventilation
- Normalized labelling provided with the equipment
- Especially designed to store acids and bases



## TECHNICAL FEATURES

Model	<b>APC112+</b>	Door closure	<b>manual closing</b>
Cabinet materials	<b>Melamine wood structure</b>	Shipping dimensions HxWxD (mm)	<b>1300 x 1200 x 800</b>
Doors color	<b>white RAL 9010</b>	Ventilation	<b>Both inlet and outlet</b>
Number of doors	<b>2</b>	Openings diameter (mm)	<b>100</b>
Doors type	<b>sliding doors</b>	Number of compartments	<b>2</b>
Hinges	<b>right &amp; left</b>	HS Code	<b>94036000</b>

## ADDITIONAL ACCESSORIES RECOMMENDED – FILTRATION OR VENTILATION BOXES

- Connecting kit (Réf : KRC)

## OPTIONS

- Containment sump built in polypropylene (Réf : E8A2)

Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Storage capacity (L)	Containment sump	Shelves	Load capacity	Weight (kg)
APC112+	1100 x 1135 x 615		100		0	20	110.00



**TRIONYX**  
DEPUIS 1971