

1- door underbench melamine safety cabinet for acids and bases

COMPLIANCE

ADVANTAGES

- Made in France
- Complies with risk prevention and safety recommendations
- Storage of corrosives

BUT ALSO...

- Adjustable leveling feet for enhance stability
- Adjustable leveling feet for enhance stability
- PP containment sumps on sliding rails
- Complies with risk prevention and safety recommendations
- Built in melamine with a great resistance to corrosive vapors
- Holes for the connection to the ventilation
- Normalized labelling provided with the equipment
- Especially designed to store acids and bases



TECHNICAL FEATURES

Model	APC32+	Door closure	manual closing
Cabinet materials	Melamine wood structure	Shipping dimensions HxWxD (mm)	870 x 600 x 800
Doors color	white RAL 9010	Ventilation	Both inlet and outlet
Number of doors	1	Openings diameter (mm)	100
Doors type	swing doors	Number of compartments	1
Hinges	right	HS Code	94036000

ADDITIONAL ACCESSORIES RECOMMENDED – FILTRATION OR VENTILATION BOXES

- Ventilation box for a connection to the outside (Réf : CDV-A)
- Connecting kit (Réf : KRC)
- Vertical recycling air box (delivered without filter) (Réf : H40+)

OPTIONS

- Containment sump built in polypropylene (Réf : E8A2)

Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Storage capacity (L)	Containment sump	Shelves	Load capacity	Weight (kg)
APC32+	630 x 635 x 578		40		0	20	30.00