

## Ventilation box for a connection to the outside

### COMPLIANCE

ISO 3744, CE

### ADVANTAGES

- CE certified
- ISO 3744 compliant
- Storage of flammables, corrosives and toxics

### BUT ALSO...

- Simple and aesthetic fitting
- Allow the extraction of the chemicals vapors contained in the cabinet
- Fan in compliance with CE specifications
- Suitable on any safety cabinet
- Easy and aesthetic installation
- Lighting control allowing to survey the proper working of the fan



### TECHNICAL FEATURES

Model	<b>CDV-A</b>	Connecting sleeve	<b>2</b>
Compliance	<b>ISO 3744, CE</b>	Electrical power (W)	<b>63</b>
Cabinet materials	<b>12/10th epoxy painted steel structure</b>	Sound power of the fan dB(A)	<b>57.00</b>

### COMPLEMENTARY PRODUCTS

- Clamp collar - diameter 100 mm (Réf : CDS100)
- Connection sleeve - diameter 100 mm (Réf : CHJ)
- Flexible sleeve in PVC - 1 meter (Réf : KL100)
- Connecting kit (Réf : KRC)

Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Weight (kg)
CDV-A	235 x 235 x 250		8.00