

Vacuum tubing transparent DN 8 mm

www.imlab.eu



Product description

The transparent vacuum tubing offers maximum performance for demanding laboratory applications. Thanks to its transparency, it enables reliable visual control of the flow, while its excellent chemical resistance guarantees a wide range of applications. The smooth inner wall effectively prevents deposits and simplifies cleaning processes. The optimum balance of flexibility and stability enables both easy installation and a secure fit at connections - even in complex set-ups. Particularly noteworthy is the extremely long service life thanks to the excellent abrasion and bending resistance. With its reinforced construction, the hose offers maximum operational reliability for vacuum applications.

Price per meter; please specify your desired length simply by selecting the quantity.

Please note: The DN 8 hose size is not suitable for applications on drying ovens - use the DN 10 size instead to achieve optimum results.

Performance features

- Price per meter, please choose desired length via quantity selection
- Excellent chemical resistance
- Smooth inner wall prevents deposits and makes cleaning easier
- Long service life due to high abrasion and bending resistance
- Maximum operational safety for vacuum applications

Ordering data

Product name	Nominal mains voltage mains frequency	Plug type	Articleno.
DN 8 mm			20686061

Items supplied

Vacuum tubing transparent DN 8 mm (sale per meter).

Technical data

Dimension B in mm	15.9mm
Dimension C in mm	7.9mm

Dimension sheet

