


KERN 474-926

Interfacekabel RS-232, Kabellengte 1,5 m



Categorie

Merk	KERN
Productcategorie	Weegschalen
Productgroep	Kabel

Ontwerp

Kabellengte	1,5 m
-------------	-------

Verpakking & verzending

Afreesbaarheid kracht [d] (N)	1 d
Verzendmethode	Pakketdienst
Afmetingen verpakking (B×D×H)	100×170×30 mm
Nettogewicht ca.	0,10 kg
Brutogewicht ca.	0,30 kg
Verzendgewicht	0,253 kg

Piktogramme