

KERN OZL 462

Tubus-trinoculaire

KERN



Zoom [Max]	4,5 x
Blikveld [Max]	4,4 mm
Monsterhoogte [Max]	140 mm

Verpakking & verzending

Afreesbaarheid kracht [d] (N)	1 d
Verzendmethode	Pakketdienst
Afmetingen verpakking (B×D×H)	385×370×240 mm
Nettogewicht ca.	1,2 kg
Brutogewicht ca.	2,0 kg

Piktogramme

STANDARD



Categorie

Merk	Optics
Productcategorie	Microscopen
Productgroep	Stereomicroscopen

Approval

CE-teken	✓
----------	---

Ontwerp

Optisch systeem	Greenough
Tubus type	Trinoculaire
Hoek van de tubus	45°
Tubus 360° rotatie	✓
Straalpadverdeling	100:0
Oculaire vast	✓
Dioptrische aanpassing [Min]	-5
Dioptrische aanpassing [Max]	5
Dioptrische aanpassing	beide zijden
Interpupilaire afstand [Min]	55 mm
Interpupilaire afstand [Max]	75 mm

Oculair

Oculair veldbreedte	HWF
Oculair vergroting	10 x
Oculair diameter	30 mm
Oculair gezichtsveld	20 mm
Gezichtspunt	High Eye Point
Objectieven vergroting	0,7 x - 4,5 x

Focus

Zoom [Min]	0,7 x
------------	-------